



xTIP[®]

By Biotix

xTIP[®] Pipette Tips

Compatible with Rainin LTS Style Pipettes

xTIP[®] Pipette Tips feature innovation beyond ergonomics. FlexFit[®], X-Resin[®], and Blade[®] technologies, developed by Biotix, are found exclusively in Biotix and xTIP brand pipette tips. Biotix technologies minimize inaccuracies associated with the number one issue in pipetting: user error. These proven features are standard on all xTIP Low Retention Pipette Tips, maximizing accuracy and precision.

*No worrying, no guesswork.
Just the best tip on the market.*

Biotix[®]

www.biotix.com



xTIP® Features

INNOVATION	BENEFITS
FLEXFIT®	FlexFit's flexible proximal end ensures confidence in a secure seal with lower insertion and ejection forces
X-RESIN®	Naturally low retention minimizes sample loss during pipetting
BLADE®	Minimize surface area on distal end of tips leading to better precision

Additional Benefits to the xTIP® Include:

- Modern tip design and high quality manufacturing, with comfort and accuracy to maximize your results.
- Alternate supplier of tips so your research continues smoothly.

Biotix® xTIP Pipette Tips, Filter Tips, Low Retention			
Volume	20 µL	200 µL	1000 µL
Sterile	Yes	Yes	Yes
Racked	R-0020-9FC	R-0200-9FC	R-1000-9FC

Biotix® xTIP Pipette Tips, Non Filter Tips, Low Retention						
Volume	20 µL		250 µL		1000 µL	
Sterile	No	Yes	No	Yes	No	Yes
Racked	R-0020-9NC	R-0020-9SC	R-0250-9NC	R-0250-9SC	R-1000-9NC	R-1000-9SC
CleanPak Reload	R-0020-9TN	R-0020-9TS	R-0250-9TN	R-0250-9TS	R-1000-9TN	R-1000-9TS
Tip Eject Reload		R-0020-9RC		R-0250-9RC		R-1000-9RC



www.biotix.com