





Biotix[®]
Fluid Innovation

xPIPETTE[™]

Biotix xTIP4 and Rainin LTS Compatible

Countless hours are spent pipetting in the lab. What impact do those long hours have on your hand, wrist and shoulders? Repetitive motions, like pipetting, can have long-lasting effects on your joints from increased strain, which may lead to repetitive stress injuries. Biotix xPIPETTE minimizes these forces to improve ergonomics in your lab.

FlexFit™ + StarStop™ = Perfect Ergonomics



Light Spring Force

- Improves ergonomics with every draw by reducing plunger forces

Secure Volume Lock

- Locked in accuracy for more reproducible results

Large Cushioned Ejector

- Silicone cushion reduces dynamic stress, lowering ejection force
- Large button design reduces mechanical stresses when you eject

Finger Hook

- Allows you to rest your hand in between pipetting cycles improving ergonomics



xPIPETTE single channel pipettes Compatible with xTIP4™ and Rainin® LTS™ pipette tips

CAT. Number	Volume Range (μL)	Volume (μL)	Accuracy (%)	Accuracy (μL)	Precision (%)	Precision (μL)
63300160 ●	0.2-2	0.2	12	0.024	6	0.012
		1	2.7	0.027	1.3	0.013
		2	1.5	0.03	0.7	0.014
63300161 ●	1-10	1	2.5	0.025	1.2	0.012
		5	1.5	0.075	0.6	0.03
		10	1	0.1	0.4	0.04
63300162 ●	2-10	2	7.5	0.15	2	0.04
		10	1.5	0.15	0.5	0.05
		20	1	0.2	0.3	0.06
63300163 ●	10-100	10	3.5	0.35	1	0.1
		50	0.8	0.4	0.24	0.12
		100	0.8	0.8	0.15	0.15
63300164 ●	20-200	20	2.5	0.5	1	0.2
		100	0.8	0.8	0.25	0.25
		200	0.8	1.6	0.15	0.3
63300165 ●	100-1000	100	3	3	0.6	0.6
		500	0.8	4	0.2	1
		1000	0.8	8	0.15	1.5