

Biotix® xPIPETTE™ and Rainin® LTS™ Compatible Pipette Tips

Perfect Ergonomic System

FlexFit®




Flexible Proximal End

- Requires less insertion force than competitors' to make a secure seal, reducing the risk of repetitive stress injury
- Tips don't fall off multichannel pipettes
- Completes the ergonomic system on LTS pipettes

StarStop™

Positive Stop Feature

- Gives consistently low ejection forces on LTS pipettes
- Aligns the tips for even liquid levels on multichannels



Clean Delivery

X-Resin®

Low Retention Surface

- Maximizes sample recovery
- Improves accuracy across your experiments

Blade®

Minimizes Surface Tension

- Eliminates hanging droplet formation
- Saves time by reducing the need for tip touch off
- Increases accuracy and precision

Biotix xTIP4

Volume μ L	Tips/Pack	Pre-Sterilized	Racked, Filtered	Racked	CleanPak Reload	Tip Eject Reload
20	960	Yes	63300931	63300005	63300010	63300015
		No		63300020	63300025	
200	960	Yes	63300001			
		No				
250	960	Yes		63300006	63300011	63300016
		No		63300021	63300026	
300	960	Yes	63300002	63300007	63300012	
		No		63300022	63300027	
1000	768	Yes	63300003	63300008	63300013	63300018
		No		63300023	63300028	
1200	768	Yes	63300004	63300009	63300014	63300019
		No		63300024	63300029	

Biotix® xPIPETTE™ and Rainin® LTS™ Compatible Pipette Tips

Structural Support (A)


- More rigid lid design for repeatable autoclave

Locking Clasp (B)

- Secures tip during transport
- Clicks to ensure rack is securely closed

Snap Card Tabs (C)

- Simply squeeze and pull up for easy one handed snap card removal with reloads



Double Ball Hinge(D)

- Easy to remove and reattach
- Opens past 180° for optimal multichannel pipetting

Structural Support (E)

- Rigid base design for repeatable autoclaving
- Stable base profile inhibits tipping of the rack when pipetting the corners

Modified Hole Design (F)

- Prevents rack lift up when removing tips at an angle



Tip Eject Reload

Convenience

- Clear casing and colored snap cards for easy tip volume identification
- Snap cards click into place
- Easily operated with one hand without contact with tips or snap card

Sustainability

- 100% recyclable packaging
- Empty racks can be autoclaved for multiple uses

Contamination Prevention

- Tips are shielded from the environment and user handling
- Rigid plastic case provides a cleaner reload than cardboard alternatives
- All Tip Eject reloads come pre-sterilized and ready to use