



Biotix[®]

xTIP4 Brochure

The Benefits of xTIP4

Experience four patented features in the lab space

Perfect Ergonomic System

FlexFit®

Flexible Proximal End

- Requires less insertion force than competitors' to make a secure seal, reducing the risk of repetitive stress injury
- Tips don't fall off multichannel pipettes
- Completes the ergonomic system on LTS pipettes

StarStop™

Positive Stop Feature

- Gives consistently low ejection forces on LTS pipettes
- Aligns the tips for even liquid levels on multichannels



BIOTIX

FlexFit®

Lower insertion force to make a secure seal

StarStop™

Eliminates high ejection forces



Clean Delivery

X-Resin®

Low Retention Surface

- Maximizes sample recovery
- Improves accuracy across your experiments

Blade®


Minimizes Surface Tension

- Eliminates hanging droplet formation
- Saves time by reducing the need for tip touch off
- Increases accuracy and precision

Biotix Tip Eject Reloads

This Tip Eject reload system offers convenience and quality while reducing your lab's waste by 30%

Structural Support (A)	Locking Clasp (B)	Snap Card Tabs (C)
<ul style="list-style-type: none">• More rigid lid design for repeatable autoclave	<ul style="list-style-type: none">• Secures tip during transport• Clicks to ensure rack is securely closed	<ul style="list-style-type: none">• Simply squeeze and pull up for easy one handed snap card removal with reloads



Double Ball Hinge (D)	Structural Support (E)	Modified Hole Design (F)
<ul style="list-style-type: none">• Easy to remove and reattach• Opens past 180° for optimal multichannel pipetting	<ul style="list-style-type: none">• Rigid base design for repeatable autoclaving• Stable base profile inhibits tipping of the rack when pipetting the corners	<ul style="list-style-type: none">• Prevents rack lift up when removing tips at an angle



Tip Eject Reload		
Convenience	Sustainability	Contamination Prevention
<ul style="list-style-type: none">• Clear casing and colored snap cards for easy tip volume identification• Snap cards click into place• Easily operated with one hand without contact with tips or snap card	<ul style="list-style-type: none">• 100% recyclable packaging• Empty racks can be autoclaved for multiple uses	<ul style="list-style-type: none">• Tips are shielded from the environment and user handling• Rigid plastic case provides a cleaner reload than cardboard alternatives• All Tip Eject reloads come pre-sterilized and ready to use

Biotix xTIP4 SKUs

Our xTIP4 tips variety in options reduce your lab's waste by 50% while providing tips with extra protection during storage

Volume μ L	Tips/Pack	Pre-Sterilized	Racked, Filtered	Racked	Tip Eject Reload
20	960	Yes	63300931	63300005	63300015
		No		63300020	
200	960	Yes	63300001		
		No			
250	960	Yes		63300006	63300016
		No		63300021	
300	960	Yes	63300002	63300007	
		No		63300022	
1000	768	Yes	63300003	63300008	63300018
		No		63300023	
1200	768	Yes	63300004	63300009	63300019
		No		63300024	



[BIOTIX.com](https://www.biotix.com)

©2021, Biotix, Inc. All rights reserved.

Biotix, xPIPETTE, xTIP4, FlexFit, StarStop, X-Resin, Blade, Tip Eject, and CleanPak are trademarks of Biotix inc.

All other brands and names contained herein are the property of their respective owners as found at [biotix.com/](https://www.biotix.com/)