

OPTIMUS RACK

By Biotix

STRUCTURAL SUPPORT

- More rigid lid design for repeatable autoclave

LOCKING CLASP

- Secures tip during transport
- Clicks to ensure rack is securely closed

FUNCTIONAL TABS

- Simply squeeze and pull up for easy one handed snap card removal with reloads

DOUBLE BALL HINGE

- Easy to remove and reattach
- Opens past 180° for optimal multi-channel pipetting

MODIFIED HOLE DESIGN

- Prevents rack lift up when removing tips

STRUCTURAL SUPPORT

- Withstands repeatable autoclaving
- Stable base profile inhibits tipping of the rack when pipetting the corners