


# Biotix uTIP Pipette Tips

Compatible with universal pipettes including Gilson®, Eppendorf®, Biohit®, and more!

## Clean Delivery

Biotix tips are molded with our proprietary low retention **X-Resin™** resulting in a uniform composition throughout the tip for ultra low retention. This low binding polymer ensures maximum sample recovery of proteins, DNA, and nanoparticles while ensuring accurate and precise measurements.

uTIP tips feature fine point precision with our **Blade™** distal end. By minimizing the surface area at the end of the tip, we have decreased the chance of burdensome hanging droplets sticking at the end of your dispense. Blade technology produces precise results while saving time by reducing the need for tip touch off.



## Ergonomic Fit

**FlexFit™** sealing technology makes an instant seal on universal pipettors. This patented feature consists of thick and thin ribs in the tip's seal area creating maximum flexibility. This flexibility allows it to easily conform to the shape of the pipette barrel and produce a secure seal.

### Biotix® uTIP™ Pipette Tips, Filtered, Low-Retention, 960 tips per pack

Volume	2µL	10µL	10µL XL	20µL	100µL	200µL	300µL	1000µL	1250µL
Sterilized	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Racked	63300153	63300154	63300041	63300042	63300043	63300044	63300045	63300046	63300047

### Biotix® uTIP™ Pipette Tips, Non-Filtered, Low-Retention, 960 tips per pack

Volume	10µL		10µL XL		200µL		250µL		300µL		1000µL		1250µL	
Sterilized	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes
Racked	63300055	63300048	63300056	63300049	63300057	63300050	63300058	63300051	63300059	63300052	63300060	63300053	63300061	63300054
CleanPak	63300069		63300070		63300071		63300072		63300073		63300074		63300075	
Bulk	63300076		63300077		63300078		63300079		63300080		63300081		63300082	