

Biotix CleanPak Reload

Compatible with universal tips

Reduce your lab's waste by 50% while providing extra protection during storage.



The Biotix CleanPak reload is easy to use and has a simple design that reduces the amount of packaging waste. The bag is hermetically sealed to better protect against contaminants during transport and while being stored.

- Color-coded snapcards by tip volume for easy identification.
- Resealable bag design for protected tip storage.
- Reduces plastic waste by 50% when compared to racked tips.

Pipette tip racks are fully autoclaveable at 121° C after transfer from reload.





Step 1: Open Bag

Tear the strip at the top of the bag and remove the protective cover.



Step 2: Transfer Refill

You will hear a click when the snapcard is secured.



Step 3: Replace Lid

Replace lid in the bag and keep it with unused inserts.



Step 4: Seal for later use

Volume μL	Pre-Sterilized	xTIP CleanPak Reload	xTIP CleanPak Packaging	uTIP CleanPak Reload	uTIP CleanPak Packaging
10	Yes				10 refills/pack
	No			75841-122	
10 XL	Yes				10 refills/pack
	No			75841-124	
20	Yes		10 refills/pack		10 refills/pack
	No	76266-118			
200	Yes		8 refills/pack		10 refills/pack
	No			75841-126	
250	Yes		10 refills/pack		10 refills/pack
	No	76266-120		76175-286	
300	Yes		10 refills/pack		10 refills/pack
	No	76266-122		75841-128	
1000			8 refills/pack		10 refills/pack
	No	76266-124		76175-386	
1200	Yes		8 refills/pack		10 refills/pack
	No	76266-126		75841-130	