



Biotix uTIP Pipette Tips

Compatible with universal pipettes including Gilson®, Eppendorf®, Biohit®, and more!

Clean Delivery

Biotix tips are molded with our proprietary low retention **X-Resin™** resulting in a uniform composition throughout the tip for ultra low retention. This low binding polymer ensures maximum sample recovery of proteins, DNA, and nanoparticles while ensuring accurate and precise measurements.

uTIP tips feature fine point precision with our **Blade™** distal end. By minimizing the surface area at the end of the tip, we have decreased the chance of burdensome hanging droplets sticking at the end of your dispense. Blade technology produces precise results while saving time by reducing the need for tip touch off.

Ergonomic Fit

FlexFit™ sealing technology makes an instant seal on universal pipettors. This patented feature consists of thick and thin ribs in the tip's seal area creating maximum flexibility. This flexibility allows it to easily conform to the shape of the pipette barrel and produce a secure seal.

Biotix® uTIP™ Pipette Tips, Filtered, Low-Retention, 960 tips per pack

Volume	2µL	10µL	10µL XL	20µL	100µL	200µL	300µL	1000µL	1250µL
Sterilized	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Racked	89425-662	89136-572	89136-576	89136-582	89136-580	76175-406	89136-568	76175-412	89408-550

Biotix® uTIP™ Pipette Tips, Non-Filtered, Low-Retention, 960 tips per pack

Volume	10µL		10µL XL		200µL		250µL		300µL		1000µL		1250µL	
	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes
Racked	10027-666	89136-570	10027-668	89136-574	10027-670	89136-578	76175-282	76175-408	10027-672	89136-566	76175-382	76175-414	10027-736	89408-548
CleanPak	75841-122		75841-124		75841-126		76175-286		75841-128		76175-386		75841-130	
Bulk	75982-122		75982-124		75982-126		76175-410		75982-128		76175-388		75982-130	