

Reduce your lab's waste by 50% while providing extra protection during storage.



The Biotix CleanPak reload is easy to use and has a simple design that reduces the amount of packaging waste. The bag is hermetically sealed to better protect against contaminants during transport and while being stored.

- Color-coded snapcards by tip volume for easy identification.
- Resealable bag design for protected tip storage.
- Reduces plastic waste by 50% when compared to racked tips.

Pipette tip racks are fully autoclaveable at 121° C after transfer from reload.



Learn more about our sustainability initiatives at [BIOTIX.com](http://BIOTIX.com)



### Step 1: Open Bag

Tear the strip at the top of the bag and remove the protective cover.



### Step 2: Transfer Refill

You will hear a click when the snapcard is secured.



### Step 3: Replace Lid

Replace lid in the bag and keep it with unused inserts.



### Step 4: Seal for later use

Volume $\mu\text{L}$	Pre-Sterilized	xTIP4 CleanPak Reload	xTIP4 CleanPak Packaging	uTIP CleanPak Reload	uTIP CleanPak Packaging
10	Yes				10 refills/pack
	No			12111021	
10 XL	Yes				10 refills/pack
	No			12111022	
20	Yes	12111376	10 refills/pack		10 refills/pack
	No	12111381			
200	Yes		8 refills/pack		10 refills/pack
	No			12111023	
250	Yes	12111377	10 refills/pack		10 refills/pack
	No	12111382		12111130	
300	Yes	12111378	10 refills/pack		10 refills/pack
	No	12111383		12111024	
1000	Yes	12111379	8 refills/pack		10 refills/pack
	No	12111384		12111136	
1200	Yes	12111380	8 refills/pack		10 refills/pack
	No	12111385		12111025	

Distributed by Fisher Scientific. Contact us today:

#### In the United States

Order online: [fishersci.com](http://fishersci.com)

Call customer service: 1-800-766-7000

#### In Canada

Order online: [fishersci.ca](http://fishersci.ca)

Call customer service: 1-800-234-7437

