

Introducing Pure Resin

Biotix is a sustainable manufacturer conscious of the impacts life science consumables have on the environment. We strive to leave nothing behind. With awareness of the impact additives from disposable plastic lab consumables have on the environment*, Pure Resin is made of 100% virgin polypropylene. No additives are released into the environment. Pure resin products also comes in 100% recyclable packaging to help you advance your research and sustainability efforts.




LEAVE NOTHING BEHIND®



Pure uTIPS

Enhance the performance and accuracy of your experiments and improve global health with Pure uTIPS. Pure utips are made with pure resin and patented technologies to improve accuracy and precision. Flexfit™ lowers necessary insertion and ejection force by ≈70% to make a secure seal, Blade™ lowers the chance of hanging droplet formation, and each tip is free of RNase, DNase, DNA, endotoxins, and nucleic acids. Same great tip and better for the environment!



Contact your local sales rep to order!

Name: _____
Email: _____
Phone: _____

Volume	Biotix Cat. No.	Description	Packaging	Filter	Sterile	Qty
2 uL	63305403	2 µL Pure uTIP™ Filter Sterile Pipette Tips	R	F	S	96 tips/rack 10 racks/pack 5 packs/case
10uL	63305404	10 µL Biotix Pure uTIP™ Filter Sterile Pipette Tips	R	F	S	
	63305412	10 µL Biotix Pure uTIP™ Sterile Pipette Tips	R	N	S	
	63305419	10 µL Biotix Pure uTIP™ Non-sterile Pipette Tips	R	N	N	
	63305426	10 µL Biotix Pure uTIP™ Bulk Non-sterile Pipette Tips	B	N	N	1,000 tips/bag 4 bags/case
10XL	63305405	10XL Biotix Pure uTIP™ Filter Sterile Pipette Tips	R	F	S	96 tips/rack
	63305413	10XL Biotix Pure uTIP™ Sterile Pipette Tips	R	N	S	10 racks/pack
	63305420	10XL Biotix Pure uTIP™ Non-sterile Pipette Tips	R	N	N	5 packs/case
	63305427	10XL Biotix Pure uTIP™ Bulk Non-sterile Pipette Tips	B	N	N	1,000 tips/bag 4 bags/case
20 uL	63305406	20 uL Biotix Pure uTIP™ Filter Sterile Pipette Tips	R	F	S	96 tips/rack 10 racks/pack 5 packs/case
100 uL	63305407	100 uL Biotix Pure uTIP™ Filter Sterile Pipette Tips	R	F	S	
200 uL	63305408	200 uL Biotix Pure uTIP™ Filter Sterile Pipette Tips	R	F	S	
	63305414	200 uL Biotix Pure uTIP™ Sterile Pipette Tips	R	N	S	
	63305421	200 uL Biotix Pure uTIP™ Non-sterile Pipette Tips	R	N	N	
	63305428	200 uL Biotix Pure uTIP™ Bulk Non-sterile Pipette Tips	B	N	N	1,000 tips/bag 4 bags/case
250 uL	63305415	250 uL Biotix Pure uTIP™ Sterile Pipette Tips	R	N	S	96 tips/rack
	63305422	250 uL Biotix Pure uTIP™ Non-sterile Pipette Tips	R	N	N	10 racks/pack 5 packs/case
	63305429	250 uL Biotix Pure uTIP™ Bulk Non-sterile Pipette Tips	B	N	N	1,000 tips/bag 4 bags/case
300 uL	63305409	300 uL Biotix Pure uTIP™ Filter Sterile Pipette Tips	R	F	S	96 tips/rack
	63305416	300 uL Biotix Pure uTIP™ Sterile Pipette Tips	R	N	S	10 racks/pack
	63305423	300 uL Biotix Pure uTIP™ Bulk Non-sterile Pipette Tips	R	N	N	5 packs/case
	63305430	300 uL Biotix Pure uTIP™ Bulk Non-sterile Pipette Tips	B	N	N	1,000 tips/bag 4 bags/case
1000 uL	63305410	1000 uL Biotix Pure uTIP™ Filter Sterile Pipette Tips	R	F	S	96 tips/rack
	63305417	1000 uL Biotix Pure uTIP™ Sterile Pipette Tips	R	N	S	10 racks/pack
	63305424	1000 uL Biotix Pure uTIP™ Non-sterile Pipette Tips	R	N	N	5 packs/case
	63305431	1000 uL Biotix Pure uTIP™ Bulk Non-sterile Pipette Tips	B	N	N	1,000 tips/bag 4 bags/case
1250 uL	63305411	1250 uL Biotix Pure uTIP™ Filter Sterile Pipette Tips	R	F	S	96 tips/rack
	63305418	1250 uL Biotix Pure uTIP™ Sterile Pipette Tips	R	N	S	10 racks/pack
	63305425	1250 uL Biotix Pure uTIP™ Non-sterile Pipette Tips	R	N	N	5 packs/case
	63305432	1250 uL uTIP™ Bulk Non-sterile Pipette Tips	B	N	N	1,000 tips/bag 4 bags/case

Legend

Packaging

R: Racked
B: Bulk

Filter

F: Filter
N: Non-Filter

Sterile

S: Sterile
N: Non-Sterile