


# uTIP Pipette Tips

Compatible with universal pipettes including Gilson®, Eppendorf®, Biohit®, and more!

## Clean Delivery

Biotix tips are molded with our proprietary low retention **X-Resin™** to ensure maximum sample recovery of proteins, DNA, and nanoparticles while ensuring accurate and precise measurements. In addition our **Blade™** distal end technology produces precise results while saving time by reducing the need for tip touch off.



## Ergonomic Fit

**FlexFit™** sealing technology makes an instant seal on various universal pipettors by easily conforming to the shape of the pipette barrel and produces a secure seal.

### Biotix® uTIP™ Pipette Tips

Biotix P/N	SLS#	Product Description	Volume	Sterilized	Filtered	Packaging Type	Qty
63300154	SLE5146	0.1-10µL Universal Pipette Tip	10µL	Yes	Yes	Rack	960 Tips/ PACK 5 PACKS / CASE
63300048	SLE5112			Yes	No	Rack	
63300055	SLE5124			No	No	Rack	
63300069	SLE5136			No	No	Reload	
63300041	SLE5148	0.1-10µL Universal Fit	10µL XL	Yes	Yes	Rack	
63300049	SLE5110			Yes	No	Rack	
63300056	SLE5122			No	No	Rack	
63300070	SLE5134			No	No	Reload	
63300042	SLE5150	1-20µL Universal Fit	20µL	Yes	Yes	Rack	
63300043	SLE5152	10-100µL Universal Fit	100µL	Yes	Yes	Rack	
63300044	SLE5154	20-200µL Universal Fit	200µL	Yes	Yes	Rack	
63300050	SLE5114			Yes	No	Rack	
63300057	SLE5126			No	No	Rack	
63300071	SLE5138			No	No	Reload	
63300045	SLE5156	20-300µL Universal Pipette Tip	300µL	Yes	Yes	Rack	
63300052	SLE5116			Yes	No	Rack	
63300059	SLE5128			No	No	Rack	
63300073	SLE5140			No	No	Reload	

## Biotix® uTIP™ Pipette Tips

Biotix P/N	SLS #	Product Description	Volume	Sterilized	Filtered	Packaging Type	Qty
63300046	SLE5158	100-1000µL Universal Pipette Tip	1000µL	Yes	Yes	Rack	960 TIPS/ PACK 4 PACKS / CASE
63300053	SLE5118			Yes	No	Rack	
63300060	SLE5130			No	No	Rack	
63300074	SLE5142			No	No	Reload	
63300047	SLE5160	100-1250µL Universal Pipette Tip	1250µL	Yes	Yes	Rack	
63300054	SLE5120			Yes	No	Rack	
63300061	SLE5132			No	No	Rack	
63300075	SLE5144			No	No	Reload	
63300151	SLE5164	10-300µL Empty Rack	10-300µL	No	No	Rack Reload	10 RACKS/PACK 5 PACKS/ CASE
63300152	SLE5166	1000-1250µL Empty Rack	1000-1250µL	No	No	Rack Reload	



BIOTIX.com