

uPIPETTE Single Channel Adjustable Volume Pipettes



Patented ejection technology reduces fatigue and makes the uPIPETTE a top choice to improve your lab's ergonomics. Designed with durable and lightweight components to produce a precise and accurate liquid transfer while maintaining hand comfort. Includes a universal shaft for flexibility to choose the right tips of your experiment.

Catalog No.	Description
30859288	0.2-2.0 μ L Single Channel, uPIPETTE
30859289	0.5-10 μ L Single Channel, uPIPETTE
30859290	2.0-20 μ L Single Channel, uPIPETTE
30859291	10-100 μ L Single Channel, uPIPETTE
30859292	20-200 μ L Single Channel, uPIPETTE
30859293	20-300 μ L Single Channel, uPIPETTE
30859294	100-1000 μ L Single Channel, uPIPETTE

- Calibrated to ISO 8655 standards and fully compatible with universal Tips
- Light spring forces produce a smooth and controlled aspiration with perfect balance
- Padded ejector button reduces the dynamic forces and vibrations your hand is exposed to during tip ejection
- Simple volume adjustment by turning the plunger

Biotix Universal Ergonomic System for flexible low force results



Biotix Universal Ergonomic System

- FlexFit seal technology creates a secure seal on contact with the pipette
- The thin, flexible walls of the seal area creates a very low force fit