


VWR®

# Universal Pipet Tips

 **vwr**™  
part of avantor

01. CLEAN DELIVERY

02. PERFECT FIT



 **VWR**®  
COLLECTION



Compatible with universal pipettors including Gilson®, Eppendorf®, Biohit®, and more!

### Secure Seal

- FLEXIBLE PROXIMAL END**
  - Gives a consistently secure seal with the widest range of industry-leading pipettors
  - Requires less insertion force than competitors' to make a secure seal, reducing the risk of repetitive stress injury
  - Tips don't fall off your multichannel pipettor
  - Makes every pipettor more ergonomic

### Clean Delivery

- LOW RETENTION SURFACE**
  - Maximizes sample recovery
  - Improves accuracy across your experiments
- MINIMIZES SURFACE TENSION**
  - Eliminates hanging droplet formation
  - Saves time by reducing the need for tip touch off
  - Increases accuracy and precision

#### VWR® Universal Pipet Tips, Low Retention, Filtered (960 Tips/Pack)

Volume	2 µL	10 µL	10 µL, XL	20 µL	100 µL	200 µL	300 µL	1000 µL	1250 µL
Sterilized	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Racked	76327-214	76322-528	76322-132	76322-134	76322-136	76322-150	76322-152	76322-154	76322-156

#### VWR® Universal Pipet Tips, Low Retention, Non-Filtered (960 Tips/Pack)

Volume	10 µL		10 µL, XL		200 µL		250 µL		300 µL		1000 µL		1250 µL	
Sterilized	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes
Racked	76322-140	76322-512	76322-142	76322-514	76322-144	76322-516	76322-146	76322-518	76322-148	76322-520	76322-524	76322-522	76322-526	76322-138
Reload	76323-394	76323-396	76323-388	76323-398	76323-390	76323-400	76323-392	-	76323-452	76323-402	76323-454	-	76323-456	76323-404

#### Empty Racks

Volume	10 - 300 µL	1000 - 1250 µL
Cat. No.	76323-406	76323-408



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by country, federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. Offers valid in countries listed above, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.