



# Biotix uTIP Pipette Tips

Compatible with universal pipettes including Gilson®, Eppendorf®, Biohit®, and more!

## Clean Delivery

Biotix tips are molded with our proprietary low retention **X-Resin™** resulting in a uniform composition throughout the tip for ultra low retention. This low binding polymer ensures maximum sample recovery of proteins, DNA, and nanoparticles while ensuring accurate and precise measurements.

uTIP tips feature fine point precision with our **Blade™** distal end. By minimizing the surface area at the end of the tip, we have decreased the chance of burdensome hanging droplets sticking at the end of your dispense. Blade technology produces precise results while saving time by reducing the need for tip touch off.



## Ergonomic Fit

**FlexFit™** sealing technology makes an instant seal on universal pipettors. This patented feature consists of thick and thin ribs in the tip's seal area creating maximum flexibility. This flexibility allows it to easily conform to the shape of the pipette barrel and produce a secure seal.

Biotix® uTIP™ Pipette Tips, Filtered, Low-Retention, 960 tips per pack

Volume	2µL	10µL	10µL XL	20µL	100µL	200µL	300µL	1000µL	1250µL
Sterilized	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Racked	1164L40	1177V54	1177V55	1177V56	1177V57	1162B04	1177V58	1162B05	1177V59

Biotix® uTIP™ Pipette Tips, Non-Filtered, Low-Retention, 960 tips per pack

Volume	10µL		10µL XL		200µL		250µL		300µL		1000µL		1250µL	
Sterilized	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes	No	Yes
Racked	1177V65	1177V60	1177V66	1177V61	1177V67	1177V62	1162B07	1162B06	1177V68	1177V63	1162B09	1162B08	1177V69	1177V64
CleanPak	1177V75		1177V76		1177V77		1162B11		1177V78		1162B13		1177V79	