

Biotix uTIP Tip Eject Reload

Compatible with universal pipettes including uPIPETTE™, Gilson®, Eppendorf®, Biohit®, and more!

Experience convenience and quality while reducing your lab's waste by 30%



This self-aligning reload solution allows you to easily reload your racks without coming into contact with the tips. The rigid 100% recycled plastic covering shields unused tips from contamination for later use.

- Available for universal compatible tips.
- 100% recycled packaging material
- Easy-to-use fast reload system
- Reduces plastic waste by 30% when compared to racked tips.
- Vertically stacked snapcards take up less space on your lab bench.

Pipette tip racks are fully autoclaveable at 121° C after transfer from reload.



Learn more about our sustainability initiatives at BIOTIX.com



Step 1: Remove the base

Remove the lower protective rack base.



Step 2: Align empty rack

Align and place Tip Eject over your empty Biotix rack.



Step 3: Push down

Push down until a “click” is heard, indicating tips are locked into rack.



Step 4: Replace stack over empty rack base

Repeat to refill empty racks as needed

Volume μL	Sterilized	uTIP Tip Eject Reload	uTIP Tip Eject Packaging
10	Yes	63305953	10 refills/pack
10XL	Yes	63305954	10 refills/pack
200	Yes	63305965	10 refills/pack
250	Yes	63305955	10 refills/pack
300	Yes	63305956	10 refills/pack
1000	Yes	63305957	8 refills/pack
1250	Yes	63305958	8 refills/pack

Volume μL	Empty Rack	Packaging
10-300	63300151	10 Rack/Pack
1000-1200	63300152	10 Rack/Pack